

Grain Vita Tester



Application:

To measure the viability of grain (germination capacity). The Grain Vita Tester works fully automatically and is microprocessor controlled.

With prepared cereal samples, the measuring instrument can be used to determine the germination capacity of a variety of different cereal species in just a few minutes with the help of a reagent fluid. The Aerne Analytic Grain Cutter is used for cutting grain grains and together with the Grain Vita Tester forms an unbeatable team.

Technical Data:

Casing: PVC white

Dimensions/mm: Depth: 280mm Width: 210 mm Height (without cover): 225 mm

Power supply: 230VAC

Scope of supply: Grain Vita Tester, 1 measuring cylinder PP 250 ml, 1 l test tube liquid (3.5 g triphenyle tetrazolium chloride) dissolved in distilled water. 1 sample holder with handle, felt and sieve. 2 plexi glass cuvettes to fill the grain halves.

Order Details:

Part:

Grain Vita Test


Order-No.

GVT

Cuvette

GVTK

Grain Cutter



Application:

Slicer for cutting Grain.

Technical Data:

Material: Aluminium, Brass
Scope of supply: Grain Cutter with insert Plate, 2 Collecting Trays

Order Details:

Part:	Order-No.
Grain Cutter	KSN
Bracket blade	KSNEK

You could also be interested in

[Mashbath](#)

