

MIRIS DMA (Dairy Analyzer)



Application:

Analysis of Dairy. Measuring Fat, Protein and Carbohydrates with MIR-Spectroscopic.

Technical Data:

Weight: 4.2 kg
Dimensions/mm: 11.5 x 26 x 37 cm
Value: 2 Digis
Power Supply: 15V, 4A DC
Duration Analysis: 60 Seconds

Order Details:

Part:	Order-No.
MIRIS DMA (Syringe)	10-02-001-DMA
MIRIS DMA (Pump)	10-01-001-DMA

MIRIS Sonicator (Homogenisator)



Application:

Homogenisator for Dairy probes before Analysis with MIRIS DMA.

Technical Data:

Maße:

52.5 x 36.5 x 24.5 cm

Dimensions:

11 kg

Input voltage 100-240 V ~ 50/60 Hz, 2.3A.

Power Supply::

Power supply Instrument

Output voltage 18 V, 100VA

Order Details:

Part:

Order-No.

SONICATOR