

Dirt Sampler SP 500 (manual)



Application:

The dirt sampler is used to test dirt in milk powder for 500 ml sample quantity. The device is attached to standard laboratory tables via a table clamp and can be adjusted in height.

The Aerne Analytic dirt sampler SP 500 is equipped with a manual bellows and a silicone drain hose.

Technical Data:

Casing: PVC perforated filter papers
Mounting: Support rod 18/8 Clamp

Order Details:

Part:	Order-No.
Dirt Sampler SP500	SP 500

Dirt Sampler SPP 500 (automatically)



Application:

For dirt testing in milk powder for 500 ml sample volume with pump and foot switch. The automatic dirt sampler is equipped with a flow pump operated by a foot switch and a silicone drain hose.

Technical Data:

Casing:	PVC perforated filter plates overlay
Power:	150 W
Dimensions:	200x 300x 500cm (LxWxH)

Order Details:


Part:

Dirt Sampler SP500 with Pump and foot switch

Order-No.

SPP 500

You could also be interested in



Racks



Waterbath Eurotherm II