

# Syringe Warmer

---



## Application:

The Aerne Analytic Syringe Warmer is applied to warm up syringes filled with baby food, medicine, etc. to a preset temperature quickly. The device is optional available with a timer. The set temperature can be programmed to values up to 75°C. The Syringe Warmer is available for a multitude of syringe types and sizes. Through the constant warming and the pre-installed adjustments the handling is easy and the daily routine in laboratories and medical practices is a relief. Syringes can be tempered exactly and be used without another checkup. The device can also be produced customized for other types of syringes or temperatures.

## Technical Data:

<b>Ambient Temperature:</b>	+5°C ... +40 °C
<b>Temperature:</b>	37°C (other target temperatures programmable) Maximum relative humidity 80 % for temperature up to 31°C, linear decreasing down to 50 % relative humidity at temperature of 40°C
<b>Humidity:</b>	
<b>Qty. of Syringe:</b>	6
<b>Syringe Type:</b>	Braun 5ml/ Luer solo
<b>Protection class:</b>	SK I
<b>Protection System:</b>	IP 40

**Performance:** 30 Watt

**Power Supply:** 230V

**Order Details:**

**Part**

Syringe Warmer for 6 syringes of Type Braun 5ml/ Luer solo

**Order-No.**

6B-ENT-230D



*“The Aerne Analytic syringe warmer relieves the daily doctor’s office a lot. The handling is very easy because of its constant warming and its inevitable adjustments. Syringes can be used exactly tempered and without anew control.”*

---

---